# Scharlau

*The wise choice*

## Aquagent®

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### Volumetric Karl Fischer titration with Scharlau Aquagent® reagents

Test measurements using Aquagent® Complet 5 and Methanol Fast

This Application Note summarizes a series of test measurements performed with an OMNIS® Titrator and Karl Fischer reagents Aquagent® Complet 5 and Methanol Fast from Scharlau.

Three series of these determinations using various water standards were carried out. The results obtained using different water standards were found to lie in a similar range. The reproducibility of the results was determined to be very good.

Using an OMNIS titration system from Metrohm and the Scharlau Karl Fischer reagents, user determinations can be carried out quickly without any decline in the reproducibility of results.

www.metrohm.com

**Titration Application Note**

### Comparison of Water Standards for The Volumetric Karl Fischer Titration

Karl Fischer (KF) titration is the method of choice for the determination of the water content in a vast variety of samples such as pharmaceuticals, petrochemical products, plastics, foods, and beverages. Compared to other analytical methods, KF titration is a simple, quick and un-expensive technique to selectively determine the water amount present in a sample.

Appropriate instrument qualification, calibration and maintenance procedures ensure correct measurement results. The qualification procedure of the titration instrument guarantees continuous the accuracy, precision, and uptime in their daily workflow. In particular, the verification of the reagent accuracy and precision of the KF titration using certified water standards is the necessary step to operate the instrument effectively.

In this application, this step is performed for the volumetric KF titration using two commonly used water standards, i.e., the 1% liquid water standard and the solid water standard of such as a bottle of liquid (1000% water content). Both standards are tested in two different KF solvents in combination with the 5 mg/mL, one-component KF titrant.

Figure 1: The computer software KF Titrate 2020 is fully controlled by the Lab Laboratory Software.

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AQ0058 Aquagent® Coulometric AG