


Hydrogen peroxide, solution 30% w/w (110 vol), EssentQ®

Identification

H₂O₂
M = 34,01 g/mol
CAS [7722-84-1]
EC number: 231-765-0
Taric code: 2847 00 00



Synonyms

Hydrogen dioxide, Hydroperoxide

Applications

analytical chemistry, oxidizing agent, bleaching agent, for pharmaceuticals synthesizing.

Specifications

assay (permanganometric).....	approx. 30 %	arsenic (As).....	max. 0,5 ppm
acidity (as H ₂ SO ₄).....	max. 0,025 %	copper (Cu).....	max. 0,001 %
chlorides (Cl).....	max. 0,001 %	lead (Pb).....	max. 0,001 %
nitrates (NO ₃).....	max. 0,001 %	iron (Fe).....	max. 5 ppm
phosphates (as PO ₄).....	max. 0,005 %	nickel (Ni).....	max. 0,001 %
sulfates (SO ₄).....	max. 0,001 %	residue on evaporation.....	max. 0,05 %

Physical data

- Density: 1,11 g/cm³
- Solub. in water: (20 °C): miscible
- Melting point: -26 °C
- Boiling point: 107 °C
- Vapour pressure: (20 °C) ~ 18 hPa
- pH(20 °C) 2 - 4

Safety - GHS

Signal Word: Danger

Hazard Statements:

H318: Causes serious eye damage.
H302: Harmful if swallowed.



Precautionary Statements:

P280: Wear protective gloves / protective clothing / eye protection / face protection.
P264: Wash thoroughly after handling.
P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310: Immediately call a POISON CENTER or doctor / physician.
P301+P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.
P501a: Dispose of contents / container in accordance with local / regional / national / international regulations.

Transport/storage

- ADR: 5.1 OC1 II • UN 2014 • HYDROGEN PEROXIDE, AQUEOUS SOLUTION
- IMDG: 5.1 II • UN 2014 • HYDROGEN PEROXIDE, AQUEOUS SOLUTION
- IATA/ICAO: 5.1 II • UN 2014 • HYDROGEN PEROXIDE, AQUEOUS SOLUTION
- PAX: 550
- CAO: 554
- 10°C - 30°C